

Download And Convert Pre-Trained Weights

Date	9/8/2022
Team ID	PNT2022TMID28486
Project Name	AI-based localization and classification of skin disease with erythema
Maximum Marks	Marks

- Using the training images located in
[`yolostructure/Data/Source_Images/Training_Images`]([/Data/Source_Images/Training_Images]) and the annotation file [`.data_train.txt`]([/Data/Source_Images/Training_Images/vottcsv-export]) which we have created in the [previous step]([/1_Image_Annotation/]) we are now ready to train our YOLOv3 detector.
- Before getting started download the pre-trained YOLOv3 weights and convert them to the Keras format. To run both steps run the download and conversion script from within the [`.yolo_structure/2_Training`]([/2_Training/]) directory
- `python Download_and_Convert_YOLO_weights.py`

